

Session Objectives


By the end of this training, you will be able to:

- ▶ Describe the different sample prep workflows
- ▶ Discuss sample prep best practices

Illumina Sequencing Workflow

1

Library Preparation



2

Cluster Generation



cBot
MiSeq

3


Sequencing



HiSeq
HiScan SQ
GA IIx
MiSeq

4

Data Analysis



ICS/RTA
CASAVA
MSR
BaseSpace

Library Preparation Overview



Nucleic Acid




Modify to proper insert size



Add adapters with sites for:

- Flow cell binding
- Sequencing primer binding


No matter the input, all libraries end up looking similar



Dual Index Library shown

The aim of the sample prep step is to obtain nucleic acid fragments with adapters attached on both ends

Selecting the Right Kit



Not sure which kit to use?

Try the “Which Illumina Sample Prep Kit is right for me” document

- TruSeq® Custom Amplicon Kit
- TruSeq Small RNA Sample Prep Kit
- TruSeq DNA Sample Prep Kit
- TruSeq Amplicon - Cancer Panel
- Nextera® DNA Sample Prep Kit
- TruSeq RNA Sample Prep Kit
- TruSeq Custom Enrichment Kit

Area of Interest	Application	Method	Input	Illumina Solution
Targeted Resequencing	<ul style="list-style-type: none"> Variant discovery DesignStudio™- generated panels 	Hybrid capture (thousands of targets)	Human gDNA with custom oligo panel	● ●
	<ul style="list-style-type: none"> Variant discovery, validation and screening DesignStudio-generated panels Pre-designed panels 	Amplicon sequencing (hundreds of targets)	<ul style="list-style-type: none"> Human gDNA FFPE DNA 	● ●
	<ul style="list-style-type: none"> Non-model population genetics Viral load/identity in blood serum High throughput drug screening Knockout screening FFPE or fresh frozen tissue tumor normal pairs High sample volume genotyping 	Amplicon sequencing (tens of targets)	Amplicons/PCR products derived from gDNA or cDNA	●
	<ul style="list-style-type: none"> Environmental testing Animal testing Plasmids Mutagenesis Vector/plasmid inserts Transfection/Infection checking 	Metagenomics (16S rRNA)	Microbial gDNA	● ●
Small Genome Sequencing	<ul style="list-style-type: none"> Microbial or viral genomes Non-model organisms Microbial or viral exome BAC/YAC screening Detecting recombination events 	<i>De novo</i> sequencing	gDNA	● ●
	<ul style="list-style-type: none"> SNP discovery Filling gaps in <i>de novo</i> assembly Detecting indel and recombination events 	Resequencing	gDNA	● ●



DNA Sample Prep Options

TruSeq DNA Sample Prep Kit Features

Highlights

- Includes index sequence by default
- Single read or paired-end compatibility
- Methylated-Cs for epigenetic applications

Sample input and indices

- 1 μ g DNA per sample
- Index up to 96 samples
- Parallel processing of up to 96 samples

Specific considerations




- Construct completed at the ligation step
- Contains internal quality controls
- Optimized for high quality genomic DNA

Suitable for:


- MiSeq
- GA IIx
- HiSeq

TruSeq DNA Sample Prep Workflow

DNA



TruSeq DNA: Two different sample prep options



Nextera Sample Prep Kit Features

Highlights

- Fastest sample prep, 90 min.
- Single read or paired-end compatibility

Sample input and indices

- 50 ng DNA per sample
- Index up to 96 samples
- Parallel processing of up to 96 samples


Specific considerations

- Construct completed at the PCR step
- Improved performance with GC regions
- Gel-free protocol

Suitable for:

- MiSeq
- GA IIx
- HiSeq

Nextera DNA Sample Prep



Nextera XT Sample Prep Kit Features

Highlights

- Fastest sample prep, 90 min.
- Optimized for small genomes, PCR amplicons and plasmids

Sample input and indices

- 1 ng DNA per sample
- Index up to 96 samples

Specific considerations

- No library quantification needed
- Sample normalization is included

Suitable for:

- MiSeq
- GA IIx
- HiSeq

Nextera XT for PCR Amplicons


Quick workflow for long or short range PCR



Short Range - CFTR




Long Range - BRCA1/2



Which version of Nextera should I use?

	Nextera	Nextera XT
Applications	Large/complex genomes (human, mammalian, plants, invertebrates)	Small genomes, amplicons, plasmids (bacteria, archea, viruses, PCR amplicons, plasmids)
DNA input	50 ng	1 ng
Post tag cleanup?	Zymo	None
Sample normalization	Manual (e.g. bioanalyzer, qPCR)	Bead-based normalization included
Indexing	96	96



Amplicon and Enrichment Sequencing Options

TruSeq Custom Amplicon Features

Highlights

- Create projects with DesignStudio
- Integrated normalization procedure
- Amplicon Viewer software

Sample input and indices

- 250 ng DNA per sample
- Up to 384 amplicons per sample
- Index up to 96 samples


Specific considerations

- Construct completed at the PCR step
- Contains internal quality controls
- Gel-free protocol

Suitable for:

- MiSeq only!

TruSeq Custom Amplicon Workflow



GATCATTAAAGATTACTTGATCCACTGATTCACGTAACCGTAACGAAACGTA
AGTAACACACITCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTA
CGTGCACACGTAACGACACTTCTGTAAACCTTAAGATTACTTGATCCACTG
CTTCTTAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAACGAAACGTA
CAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTA
ATAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTA
GATCATTAAAGATTACTTGATCCACTGATTCACGTAACCGTAACCGTAACCGTAACCGTA

TTAAGATTACTTGATCCACTGATTCACGTAACCGTAACCGTAACCGTAACCGTA
ACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAACCGTA
CGTGCACACGTAACGACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTA
CTTCTTAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAACCGTAACCGTA
CAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTA
ATAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTAACCGTA
GATCATTAAAGATTACTTGATCCACTGATTCACGTAACCGTAACCGTAACCGTAACCGTA

TruSeq Amplicon – Cancer Panel Features

Highlights

- Most relevant cancer loci in a single panel
- Suited for FFPE samples
- Integrated normalization procedure
- Amplicon Viewer software

Sample input and indices

- 250 ng DNA per sample
- 48 genes are targeted with 212 amplicons
- Index up to 96 samples

Specific considerations

- Construct completed at the PCR step
- Contains internal quality controls
- Gel-free protocol

Suitable for:


- MiSeq only!

TruSeq Amplicon – Cancer Panel

Pre-screen Samples

Assay Biochemistry

Sequencing and Analysis



FFPE QC Kit



TruSeq Amplicon - Cancer Panel



Automated Sequencing and Variant Calling

TruSeq Enrichment Sample Prep Kit Features

Highlights

- Targeted resequencing of human exome and associated regulatory regions
- Create Custom Enrichment with DesignStudio

Sample input and indices

- 500ng TruSeq DNA library per enrichment
- Index up to 6 samples per Exome Enrichment
- Index up to 12 samples per Custom Enrichment




Specific considerations

- Contains internal quality controls
- Gel-free protocol

Suitable for:

- GA IIx
- HiSeq
- MiSeq (for certain Custom Enrichments)

TruSeq Enrichment Workflow Overview



Repeat A through D

PCR Amplify

Ready for clustering

Targeted resequencing to focus on subset of genome

- Exome
Gene encoding sequences
- Custom
Researcher/user-defined regions

Nextera Enrichment Sample Prep Kit Features

Highlights

- Targeted resequencing of human exome and associated regulatory regions
- Create Custom Enrichment with DesignStudio

Sample input and indices

- 50 ng DNA input
- Index up to 12 samples per Enrichment

Specific considerations

- Contains internal quality controls
- Gel-free protocol

Suitable for:






- GA IIx
- HiSeq
- MiSeq (for certain Custom Enrichments)

Nextera Enrichment Workflow

Targeted resequencing to focus on subset of genome

Exome
Gene encoding sequences

Custom
Researcher/user-defined regions



TruSeq Enrichment or Nextera Enrichment?

	TruSeq Enrichment	Nextera Enrichment
Library prep input	1ug	50ng
Library prep time (hrs)	12	3
Prep & Enrich kits	Separate	Bundled
Total workflow (days)	3.5	2.5
Pre-enrich pooling (E/C)	6/12	12/12
Dual-indexing support	No (Up to 24 single-indexes)	Yes (Up to 96 dual-indexes)
% Enrichment (aligned/+padded)	55%/65%	55%/65% (60%/70% typical)
% Region Dropout (E/C)	<1%/<5%	<1%/<5%
Cov Uniformity (0.2X Mean)	≥80%	≥80%
Cov Uniformity (0.5X Mean)	≥60%	≥60%

E = exome; C = custom



RNA Sample Prep Options

TruSeq RNA Sample Prep Kit Features

Highlights

- Includes index sequence by default

Sample Input and indices

- 0.1-4ug of total RNA or 10-400ng of mRNA
- Up to 24 indices
- Parallel processing of up to 96 samples

Specific considerations


- Construct completed at the ligation step
- Contains internal quality controls
- Gel-free protocol

Suitable for:


- MiSeq
- GA IIx
- HiSeq

TruSeq RNA Workflow Overview

RNA to cDNA



TruSeq DNA and RNA Workflow Overview



TruSeq small RNA Sample Prep Kit Features

Highlights

- Includes index sequence by default

Sample Input and indices

- 1 ug (total RNA) or 10-50 ng (purified miRNA)
- Up to 48 indices
- Parallel processing of up to 96 samples


Specific considerations

- Construct completed at the PCR step

Suitable for:

- MiSeq
- GA IIx
- HiSeq

TruSeq Small RNA Workflow Overview






Sample Prep Best Practices

Quantification Methods of Sequencing Libraries


UV- spectrophotometer
Nanodrop

- Detects all nucleic acids nonspecifically
- Contaminants elevate values



Bioanalyzer 2100

- Accuracy highly dependent on dilution and sample handling
- Can narrow in on library peak



Fluorescence-based
ds-DNA assay
Qubit or PicoGreen

- Detects all double-stranded DNA
- Does not discriminate incomplete libraries



qPCR


- Specifically measures full-length libraries
- Detection very sensitive



Recommended Library QC methods

		Quality Check on Bioanalyzer 2100	Recommended Quant Method	Comments
Sample Prep Method	TruSeq DNA (WGS)	✓	qPCR	
	TruSeq DNA (enrichment)	✓	Qubit	Qubit is sufficient when used after DNA sample prep and before enrichment
	TruSeq RNA	✓	qPCR	
	TruSeq Enrichment	✓	qPCR	
	TruSeq small RNA	✓	Bioanalyzer	Bioanalyzer is sufficient because reduced precision needed due to: <ul style="list-style-type: none"> • low concentration libraries • short read lengths
	Nextera/XT DNA	✓	Qubit	
	TruSeq Custom Amplicon	No	None	Built-in Normalization step Final library is single stranded

Summary





Questions?