



Troubleshooting

© 2011 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CSPPro, DASL, Eco, Genetic Energy, GALx, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

illumina[®]



Learning Objectives

After completing this lesson, you will be able to:

- ▶ Troubleshoot initializing errors
- ▶ Troubleshoot fluidic issues
- ▶ Troubleshoot the vacuum system

Agenda

Initializing
Troubleshooting

Fluidics
Troubleshooting

Vacuum
Troubleshooting






Software initialization troubleshooting

Symptoms:

Error message: Software was unable to initialize internal hardware devices





Fluidics Troubleshooting

Symptoms:

Bubbles visible in the flow cell/Waste volume too low

Reposition the FC and confirm that the wholes are facing down



Check volume in reagent bottles



Perform an instrument wash



Check gaskets

Check sippers


Fluidic System check

Wash the system

- Collect the waste of 8 tubes per flow cell used
- Place the waste line for each flow cell into its own 15 mL tube
- Run a water wash
- Check the waste volume for each waste line (~4 mL per tube)

Perform Fluidic check

- Visually inspect the flow cell for bubbles and leaks near the manifold



- A** Flow Cell Waste Lines for Reagent Positions 1–8
- B** Condensation Pump Tubing (Do not remove)
- C** Paired-End Priming Pump Tubing (Do not remove)



Vacuum Troubleshooting

Symptoms

Flow Cell Holder Switch


- Has the wrong color in position 1 or 2

Flow Cell Holder






- Is not holding vacuum
- You are able to lift a flow cell off the stage when the vacuum switch is engaged

Poor or No Vacuum

- Affects both flow cell holders if the vacuum pump is defective



HiSeq System Flow Cell Holder Switch States

Switch Position and Color	Description
	Position 0 – off. Switch not illuminated. No vacuum to the flow cell holder.
	Position 1 – vacuum engaged. Solid orange indicates that inadequate vacuum is being provided.
	Position 1 – vacuum engaged. Blinking orange indicates that vacuum is being provided, but still inadequate.
	Position 1 – vacuum engaged. Blinking green indicates that vacuum pressure is good. Switch can be moved to third position (raise manifolds).
	Position 2 – manifolds raised. Solid green indicates good vacuum pressure.


Troubleshoot the flow cell holder area

Examine and clean the flow cell holder

- Remove flow cell and check vacuum ports and groove

Check if one flow cell has vacuum

- indicates the vacuum pump is working



A Vacuum ports

B Groove in flow cell holder through which vacuum is applied



Questions?