PyroMark™ Q24

Installation and Safety

Table of Contents

Legal	1
Intended Use – United States	1
Intended Purpose – Europe	1
Warranty and Liability	1
Limitation of Warranty and Liability	1
Trademarks and Patents	2
Other Trademarks	2
Copyright	2
Install PyroMark™Q24 System	3
Install and Move PyroMark™Q24 Instrument	3
Delivered Items	3
Site Requirements	4
Install the Instrument	4
Move the Instrument	5
Upgrade the Instrument Software	6
Install and Uninstall PyroMark™Q24 Software	7
Delivered Items	7
Computer Requirements	7
Install or Upgrade PyroMark™Q24 Software	7
Uninstall PyroMark™Q24 Software	8
Install and Move PyroMark™Q24 Vacuum Prep Workstation	9
Delivered Items	9
Site Requirements	10
Install the Vacuum Prep Workstation	10
Move the Vacuum Prep Workstation	11

Installation Test	12
Instrument – Function Test	12
Vacuum Prep Workstation – Function Test	12
User Documentation	13
Installation and Safety	13
User Manual	13
Quick Guides	13
Online Help	13
Instructions for Reagents and Consumables	14
Service Agreements	14
Safety	15
Safety Requirements	15
Symbols Used on Pyrosequencing Product Labels	17
Safety Label on the Instrument	18
Symbols Used in the User Documentation	18
WEEE Compliance Statement	18
Safety in Other Languages	20
Index	46
Contact Information	48
Sales Offices and Distributors	48
Technical Support	48
Revision History	48

Legal

Intended Use – United States

PyroMark™Q24 System is designed for Laboratory Use Only, which means it may be used for either research purposes or by high complexity CLIA certified laboratories.

All operations must be performed according to instructions provided through dialogs appearing on the screen, the manuals, instructions for use, and by Pyrosequencing AB's technical support staff, and within limits set by the technical specifications.

Intended Purpose – Europe

PyroMark™Q24 System is a system for detecting changes in specified variable positions in human DNA that may have clinical relevance.

PyroMark™Q24 System is available for research and, in certain European countries*, for *in vitro* diagnostic applications. PyroMark™Q24 System meets the requirements of Annex III of the European Directive for In Vitro Diagnostic Medical Devices 98/79/EC.

For *in vitro* diagnostic medical use, PyroMark™Q24 System may only be operated (i) by personnel who have received special education and training with regard to procedures utilizing *in vitro* diagnostic medical devices, and (ii) accredited medical testing laboratories.

All operations must be performed according to instructions provided through dialogs appearing on the screen, the manuals, instructions for use, and by Pyrosequencing AB's technical support staff, and within limits set by the technical specifications.

* For more information, see www.biotagebio.com

Warranty and Liability

Pyrosequencing AB warrants for a period of twelve (12) months from the date of shipment ("Warranty Period") that PyroMark™Q24 System (i) meets the published specifications, and (ii) is free from defects in material and workmanship under normal use and service and when used in compliance with the applicable operating instructions. Pyrosequencing AB's sole liability and Buyer's exclusive remedy for a breach of this warranty is limited to replacement, repair or refund at the sole option of Pyrosequencing AB. This warranty does not apply to any consumable items included in PyroMark™Q24 System such as, but not limited to, filter probes, tubing, fittings, o-rings and gaskets.

Limitation of Warranty and Liability

PYROSEQUENCING AB HEREBY EXPRESSLY DISCLAIMS, AND BUYER HEREBY EXPRESSLY WAIVES, ANY WARRANTY REGARDING RESULTS OBTAINED THROUGH THE USE OF THE SOFTWARE, INCLUDING WITHOUT LIMITATION ANY CLAIM OF INACCURATE, INVALID, OR INCOMPLETE RESULTS. ALL OTHER WARRANTIES, REPRESENTATIONS, TERMS AND CONDITIONS (STATUTORY, EXPRESS, IMPLIED OR OTHERWISE) AS TO QUALITY, CONDITION, DESCRIPTION, MERCHANTABILITY, FITNESS FOR PURPOSE OR NON-INFRINGEMENT (EXCEPT FOR THE IMPLIED WARRANTY OF TITLE) ARE HEREBY EXPRESSLY EXCLUDED.

PYROSEQUENCING AB SHALL IN NO EVENT BE LIABLE FOR ANY INDIRECT OR CONSEQUENTIAL, OR PUNITIVE DAMAGES OF ANY KIND FROM ANY CAUSE ARISING OUT OF THE SALE, USE OR INABILITY TO USE ANY SOFTWARE, INCLUDING WITHOUT LIMITATION, LOSS OF PROFITS, GOODWILL OR BUSINESS INTERRUPTION, EVEN IF PYROSEQUENCING AB HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

Trademarks and Patents

Biotage, 1-Point Support, Pyrosequencing, Pyrogram, PyroMark, and  are trademarks owned by Biotage AB. Pyrosequencing technology is covered by patents including patents US 5405746, US 6210891, US 6258568, US 6828100, US 7078168, EP 0932700, EP 0946752, EP 1495137, JP 3510272 and JP 3533223, and patent applications owned by Biotage AB.

In view of the risk of trademark degeneration, authors intending to use the trademark designations are respectfully requested to acknowledge the trademark status of the products at least once in each article.

Other Trademarks

Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

Copyright

Pyrosequencing AB reserves the right to make changes to the information contained herein without prior notice. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the expressed written permission of Pyrosequencing AB.

© Copyright 2008 Pyrosequencing AB. All rights reserved.

Pyrosequencing AB

Kungsgatan 76
SE-753 18 Uppsala
SWEDEN

Phone: +46 18 56 59 00

Fax: +46 18 59 19 22

E-mail: pyrosequencing@biotage.com

Website: www.biotagebio.com

Pyrosequencing AB is a wholly-owned subsidiary of Biotage AB. Biotage AB is listed on the OMX Nordic Exchange (ticker BIOT), www.biotage.com.

Install PyroMark™Q24 System

WARNING! You must observe all safety requirements on page 15 - 16 when installing and operating PyroMark™Q24 System. Failure to install or use the system in a manner specified by Pyrosequencing may result in personal injury and/or equipment damage.

Note: We recommend that the boxes and packing materials are kept by the customer in case an item needs to be returned for service or moved to another location.

Install and Move PyroMark™Q24 Instrument

WARNING! The total weight of the package including the instrument is 35 kg. Use suitable lifting equipment when moving the package.

WARNING! The instrument weighs 28 kg. Two persons are required when lifting the instrument.

WARNING! Use only power supplies and cords supplied with your instrument. The mains plugs shall be easily accessible in case of the need of quickly disconnecting the instrument from the mains power.

Delivered Items

Check the delivery for completeness and for damage during transportation. The following items are delivered:

- PyroMark™Q24 Instrument
- PyroMark™Q24 Instrument Software (pre-loaded)
- PyroMark™Q24 Cooling Device (pre-installed)
- 12 VDC power supply
- 24 VDC power supply
- Power cord (2)
- USB memory stick (2)
- Installation and Safety booklet (this publication)
- Quick Guide
- Declaration of Conformity document
- Test report
- Packing list
- PyroMark™Q24 Software package (see page 7)
- PyroMark™Control Oligo (shipped in a separate package)

Note: PyroMark™Q24 User Manual is available on the PyroMark™Q24 Software CD.

Site Requirements

Before installing the instrument, the installation site should be prepared as follows:

Bench Flat, stable, and able to support the weight of the instrument (28 kg).

At least 70 x 60 cm (W x D) free area is required. The free space above the bench should be at least 70 cm to accommodate the instrument with the lid open, and for ventilation. Allow at least 10 cm of free space behind the instrument for cabling.

Environment Temperature: 15°C to 32°C*

Humidity: 20% to 90% relative humidity

Draft-free location not close to window.

Keep the instrument out of direct sunlight.

* We recommend that the instrument is run in normal room temperature (20°C to 25°C) for optimal performance

Input voltage and current 100–240 VAC, 47–63 Hz, 1.1–0.45 A

The instrument must be connected to protective earth.

Note: In the event of a power failure, data obtained will still be stored on the instrument's internal memory for later retrieval.

Install the Instrument

1. Before installing the instrument, the installation site should be prepared according to the site requirements above.
2. Place the instrument on a bench (dimensioned according to the site requirements and that are able to support the weight, 28 kg). Lift the instrument by the base plate. Two persons are needed.
3. Remove the transport strap from under the instrument lid (see image below) using a pair of scissors.
4. Ensure that the power switch at the rear of the instrument is in position 0 (off); see image below.
5. Connect the two power supplies to the instrument. Place the supplies in a ventilated area, on the bench or on the floor. Do not cover the supplies with papers or other items that will reduce cooling.

See page 19 for larger image.

6. Connect the two power supplies to properly grounded (earthed) mains outlets with the correct mains voltage and frequency using the supplied power cords. The mains plugs shall be easily accessible in case of the need of quickly disconnecting the instrument from the mains power.
7. At the rear of the instrument, press the light button and check that the coolant level is visible in the window (see image). If not, please contact 1-Point Support before proceeding (see contact details on page 48).
8. Set the correct date and time in the instrument software; see the *System Administration* section in PyroMark™Q24 User Manual.
9. Ensure that the instrument is working properly by using PyroMark™Control Oligo (one package is supplied with the system). Set up the run according to the instructions supplied with the product.

For instructions on how to use the instrument, please read PyroMark™Q24 User Manual (see the *User Documentation* section on page 13).

Move the Instrument

Should you need to move the instrument within the laboratory or to a new site, follow the instructions below.

1. When the instrument is not processing, select **Shutdown** in the instrument's main menu (using the ▲ and ▼ screen buttons) and press **Ok**.

2. When the message *It is now safe to turn off the instrument* appears, turn off (O) the instrument. The power switch is located at the rear of the instrument.
3. Disconnect the instrument from the mains power. There are two mains plugs.
4. Disconnect the two power supplies from the instrument.
5. To move the instrument within the laboratory, carefully place the instrument on a trolley and move to the new room. Lift the instrument by the base plate. Two persons are needed.

6. To transport the instrument to another site:
 - a. Clean the instrument according to the instructions in the *Maintenance* section in PyroMark™Q24 User Manual.
 - b. Secure the dispensation unit to the metal bar (behind the dispensation unit) using a transport strap; see images below.
 - c. Pack the instrument properly, preferably by using the box and packing materials from the original delivery of the instrument. Lift the instrument by the base plate. Two persons are needed.
 - d. Use suitable lifting equipment when moving the package.
7. Install the instrument at the new location according to the instructions on page 4.

Upgrade the Instrument Software

If you have received a software upgrade from Pyrosequencing, upgrade your instrument as described below.

1. Save the upgrade files on one of the USB memory sticks supplied with the system. The files should be saved in a folder called *Upgrade* at the root of the memory stick.
2. When the instrument is not processing, plug the USB memory stick into the USB port at the front of the instrument. Leave it mounted until the upgrade is finished.
3. Select **Administration** in the main menu (using the ▲ and ▼ screen buttons) and press **Ok**.
4. Select **Upgrade Software** and press **Ok**.
5. Follow the instructions on the screen.

Install and Uninstall PyroMark™Q24 Software

Delivered Items

Check the delivery for completeness and for damage during transportation. The following items are delivered:

- PyroMark™Q24 Software CD
- Quick Guide

The PyroMark™Q24 Software package is normally included in the instrument package.

Computer Requirements

The office computer used to set up runs and analyze data should be a personal computer with the following recommended minimum specifications:

Operating system	Microsoft Windows XP or Windows Vista, English versions
Processor	Intel Pentium 4, 3 GHz
RAM	1 GB
Hard drive	100 MB free space
Graphics card	Supporting the resolution of the monitor
Monitor	1280 x 1024 pixels
Pointer device	Mouse or similar
Interfaces	USB-port and CD-ROM

Install or Upgrade PyroMark™Q24 Software

Note: If your computer is connected to a network, network policy settings might prevent you from completing this procedure. For more information, contact your system administrator.

1. Ensure that the computer meets the minimum requirements above and that you are logged on as an administrator or a member of the Administrators group on your computer.
2. Close any programs running on the computer.
3. Insert the PyroMark™Q24 Software CD into the CD-drive.

4. In the CD menu that automatically appears, click **Install PyroMark Q24 Software**. If the CD menu does not appear:
 - a. Select **(My) Computer** in the Windows **Start** menu.
 - b. Right-click the CD-drive with the PyroMark™Q24 Software CD and select **Explore**.
 - c. Double-click the file *autorun.exe*.
5. Follow the instructions that appear in the Setup Wizard. When the software has been successfully installed, you will return to the CD menu.
6. To view some of the analysis reports, you must have Adobe Acrobat Reader installed on your computer. To install Adobe Acrobat Reader on your computer, click **Install Acrobat Reader** in the CD menu and follow the instructions that appear in the Setup Wizard. You can also download Adobe Acrobat Reader at www.adobe.com.
7. When you are done installing, click **Exit Setup** in the CD menu.
8. Please use Windows Update (www.update.microsoft.com) to check for any critical updates to the .NET Framework 2.0.
9. To start the software, select **All Programs | Biotage | PyroMark Q24** in the Windows **Start** menu. For instructions on how to set up and analyze a run, see the online help (press the **F1** key when in the software).

Uninstall PyroMark™Q24 Software

1. Select **Control Panel** in the Windows **Start** menu.
2. In the **Control Panel**, click **Add or Remove Programs** (XP) or **Uninstall a program** under the Programs category (Vista).
3. In the list of programs, select **PyroMark Q24**.
4. Click **Remove** (XP) or **Uninstall** (Vista).
5. Repeat steps 3 through 4 for PyroMark Launcher.

Install and Move PyroMark™Q24 Vacuum Prep Workstation

WARNING! The total weight of the package including the vacuum prep workstation is 11 kg. Use suitable lifting equipment when moving the package.

WARNING! The mains plug shall be easily accessible in case of the need of quickly disconnecting the vacuum pump from the mains power.

Delivered Items

Check the delivery for completeness and for damage during transportation. The following items are delivered:

- PyroMark™Q24 Vacuum Prep Tool
- PyroMark™Q24 Vacuum Prep Worktable
- PyroMark™Q24 Troughs (12)
- PyroMark™Q24 Plate Holder (2)
- Pump
- Waste container and cap with tubing
- Waste filter (2)
- Tubing
- Allen key
- Installation and Safety booklet (this publication)
- Quick Guide
- Declaration of Conformity document
- Test report
- Packing list

Note: PyroMark™Q24 User Manual is available on the PyroMark™Q24 Software CD.


Site Requirements

Before installing the vacuum prep workstation, the installation site should be prepared as follows:

Bench	Flat, stable, and able to support the weight of the vacuum prep workstation (11 kg including all parts and with filled troughs and a full waste container). At least 35 x 70 cm free area is required. The free space above the bench should be at least 40 cm to accommodate working space, and for ventilation.
Environment	Temperature: 15°C to 32°C Humidity: 20% to 90% relative humidity Normal laboratory conditions. Adequate ventilation.
Input voltage and current	100 VAC, 50/60 Hz, 1.7/1.4 A, or 115 VAC, 60 Hz, 1.5 A, or 230 VAC, 50 Hz, 0.6 A

Install the Vacuum Prep Workstation

1. Before installing the vacuum prep workstation, the installation site should be prepared according to the site requirements above.
2. Place the workstation on a bench (dimensioned according to the site requirements and that are able to support the weight, 11 kg).
3. Cut the tubing into three pieces and connect the tubing according to the image below. Ensure that you connect the tubing to the pump's Vacuum fitting. (See page 19 for larger image.)
4. Connect the vacuum pump to a mains outlet with the correct mains voltage and frequency (check the label underneath the pump). The mains plug shall be easily accessible in case of the need of quickly disconnecting the vacuum pump from the mains power.
5. Place six of the troughs supplied with the workstation on the worktable.
6. Ensure that the filter probes are working properly by performing the function test on page 12.

For instructions on how to use the workstation, please read PyroMark™ Q24 User Manual (see the *User Documentation* section on page 13).

Move the Vacuum Prep Workstation

WARNING! Irritant. The Denaturation Solution, which is used with the workstation, contains sodium hydroxide that is irritating to eyes and skin. Empty the waste container and tubing and discard any solutions left in the troughs before moving the workstation (see instructions below). Use adequate ventilation. MSDS can be downloaded at www.biotagebio.com.

Should you need to move the vacuum prep workstation within the laboratory or to a new site, follow the instructions below.

Note: Observe all federal, state, and local environmental regulations for disposal of laboratory waste.

Required items: Impermeable gloves, chemical safety goggles, protective clothing, high purity water (Milli-Q 18.2 MΩ x cm, www.millipore.com, or equivalent), a beaker, and lint-free cloths.

1. Ensure that the vacuum switch is closed (OFF) and the vacuum pump is turned off (O).
2. Disconnect the vacuum pump from the mains power.
3. Discard any solutions left in the troughs.
4. Empty the tubing; pull off one end of the tubing and empty it of any liquid waste using an empty beaker.
5. Empty the waste container. Note that the cap can be removed without disconnecting the tubing.

6. To move the workstation within the laboratory, carefully place the workstation on a trolley and move to the new room.
7. To transport the workstation to another site:
 - a. Disconnect the tubing.
 - b. Rinse the waste container, tubing, and troughs with high purity water (Milli-Q 18.2 MΩ x cm, www.millipore.com, or equivalent). Leave them to dry.

- c. Clean the worktable and tool, except for the filter probes, using a clean, lint-free cloth moistened with water or a mild detergent. Do not touch the tips of the filter probes.
 - d. Wipe the worktable and tool, except for the filter probes, dry using a clean, lint-free cloth.
 - e. When all parts are dry, pack the workstation properly, preferably by using the box and packing materials from the original delivery of the workstation.
8. Install the workstation at the new location according to the instructions on page 10.

Installation Test

Instrument – Function Test

Before using the instrument, check that it is working properly by using PyroMark™Control Oligo (one package is supplied with the system). Set up the run according to the instructions supplied with the product. To order PyroMark™Control Oligo, please contact your local sales representative (www.biotagebio.com).

Vacuum Prep Workstation – Function Test

Before using the vacuum prep workstation, check that the filter probes are working properly by performing the function test below.

Required items: A PCR plate, a trough, and high purity water (Milli-Q 18.2 MΩ x cm, www.millipore.com, or equivalent).

1. Add 100 µl of high purity water to each well of a PCR plate.
2. Fill a trough with 70 ml of high purity water.
3. Start the vacuum pump (I).
4. Apply vacuum to the vacuum prep tool by opening the vacuum switch.
5. Lower the filter probes into the trough with high purity water and hold them there for approximately 20 seconds. Ensure that the water is transported to the waste container, i.e. that vacuum is applied. If not, check the connections.
6. Lower the filter probes into the PCR plate and check that the water is aspirated evenly in all wells and that all wells are empty after a maximum of 10 seconds.
7. If the wells are not empty after 10 seconds, repeat the procedure from step 1. If the test fails twice, replace the failing filter probes (see the *Maintenance* section in PyroMark™Q24 User Manual).

User Documentation

User documentation in other languages can be downloaded at the Pyrosequencing Technical Support website.

Installation and Safety

PyroMark™Q24 Installation and Safety contains safety requirements and instructions for installing, moving, and re-installing the system.

A printout of the document is supplied with the system. A PDF version of the document is available on the PyroMark™Q24 Software CD and on the Pyrosequencing Technical Support website, but can also be accessed by selecting **All Programs | Biotage | PyroMark Q24 Documentation | PyroMark Q24 Installation and Safety.pdf** in the Windows **Start** menu.

User Manual

PyroMark™Q24 User Manual provides a description of the system, technical specifications, instructions for operation and maintenance, and guides for assay design and validation and troubleshooting.

The manual is available as a PDF document on the PyroMark™Q24 Software CD and on the Pyrosequencing Technical Support website, but can also be accessed by selecting **All Programs | Biotage | PyroMark Q24 Documentation | PyroMark Q24 User Manual.pdf** in the Windows **Start** menu.

Quick Guides

The three quick guides contain short descriptions on how to operate PyroMark™Q24 Instrument, PyroMark™Q24 Software, and PyroMark™Q24 Vacuum Prep Workstation.

Printouts of the guides are supplied with the system (one printout per guide). A PDF version of the guides is available on the PyroMark™Q24 Software CD and on the Pyrosequencing Technical Support website, but can also be accessed by selecting **All Programs | Biotage | PyroMark Q24 Documentation | PyroMark Q24 Quick Guide.pdf** in the Windows **Start** menu.

Online Help

The online help included with PyroMark™Q24 Software contains user instructions for PyroMark™Q24 Software.

Context-sensitive help is accessed by pressing the **F1** key when in a dialog or window in PyroMark™Q24 Software. A PDF version of the online help is available on the Pyrosequencing Technical Support website, but can also be accessed by selecting **All Programs | Biotage | PyroMark Q24 Documentation | PyroMark Q24 Online Help.pdf** in the Windows **Start** menu.

Instructions for Reagents and Consumables

Detailed instructions for use are included with PyroMark™Q24 Gold Reagents, PyroMark™Q24 Cartridge, PyroMark™Control Oligo, and PyroMark™Q24 Validation Oligo. MSDS for PyroMark™Q24 Gold Reagents, PyroMark™Control Oligo, PyroMark™Q24 Validation Oligo, and sample preparation solutions can be downloaded at the Pyrosequencing Technical Support website.

Service Agreements

Service agreements are available for PyroMark™Q24 Instrument. Please contact your local sales representative or refer to Pyrosequencing website (www.biotagebio.com) for further details.

Safety

Safety Requirements

WARNING! You must observe all safety requirements when installing and operating PyroMark™Q24 System. Failure to install or use the system in a manner specified by Pyrosequencing may result in personal injury and/or equipment damage.

- PyroMark™Q24 System shall be installed in accordance with the *Install PyroMark™Q24 System* section on page 3.
- The mains plugs shall be easily accessible in case of the need of quickly disconnecting the equipment from the mains power.
- PyroMark™Q24 System shall be used in accordance with the instructions provided through dialogs appearing on the screen, the manuals, instructions for use, and by 1-Point Support, and within limits set by the technical specifications (see PyroMark™Q24 User Manual). For more information on user documentation, see page 13.
- Follow all generally accepted laboratory safety procedures and applicable laws and regulations. Use appropriate safety precautions when handling biological material.
- The Denaturation Solution, which is used with the vacuum prep workstation, contains sodium hydroxide that is irritating to eyes and skin. Use adequate ventilation. Use impermeable gloves, chemical safety goggles, and protective clothing. MSDS can be downloaded at www.biotagebio.com. Always empty the waste container and discard any sodium hydroxide left in the troughs at the end of a working day (see PyroMark™Q24 User Manual).
- Use only power supplies and cords supplied with your instrument. Power supplies and cords should be inspected periodically and replaced if damaged or altered. Note: The ground prong on the mains plugs must not be removed and the plugs should only be connected to grounded outlets as per local and national regulations. Use of power supplies that disable the grounding connections to the instrument is dangerous and prohibited.
- Always keep the instrument lid closed while the instrument is processing. An audible warning signal will alert you if the lid is opened when it is not safe.
- Do not remove the cover panels — there are no user serviceable parts inside. If you believe there is a problem with your system, please contact 1-Point Support immediately (www.biotagebio.com).
- Always shut down and disconnect the equipment from the mains power before any maintenance is performed. Note that the power switch does not disconnect the equipment completely from mains power.

- The instrument weighs 28 kg. Two persons are required when lifting the instrument.
- If the system has been damaged or does not function properly, shut down the system and contact 1-Point Support immediately (www.biotagebio.com).
- Do not touch the sharp needles at the bottom of the reagent cartridge.
- Do not touch the hot surfaces on the plate holder and the heating block.
- Do not connect or use any consumables, accessories, or external equipment other than what is specified in the user documentation. Please refer to Pyrosequencing website (www.biotagebio.com), or contact your local sales representative, for the latest information on consumables and accessories.

Symbols Used on Pyrosequencing Product Labels

Explanation of symbols used on Pyrosequencing product labels:

Symbol Referent

In accordance with all the essential requirements of all applicable European product directives; see the Declaration of Conformity document supplied with the product

In accordance with both U.S. and Canadian safety standards; see the Declaration of Conformity document supplied with the product

In accordance with FCC (Federal Communications Commission) standards; see the Declaration of Conformity document supplied with the product

Subject to the Waste Electrical and Electronic Equipment (WEEE) Directive; see page 18

In vitro diagnostic (IVD) medical device

Catalog or reorder number

Batch code or lot number

Use by (+ date)

Temperature limitation

Consult accompanying user documentation

Manufacturer

Date of manufacture

Serial number

Warning, consult accompanying user documentation

Irritant substance

Safety Label on the Instrument

The following safety label is attached to PyroMark™Q24 Instrument:

Pinch and impact warning due to moving parts.

Keep the instrument lid closed while the instrument is processing. An audible warning signal will alert you if the lid is opened when it is not safe.

Note: For information on the symbols used on the product label at the rear of the instrument, see page 17.

Symbols Used in the User Documentation

Warning or Caution:

If the symbol is accompanied by “**WARNING!**”, it means that you must follow an instruction in detail to avoid personal injury.

If the symbol is accompanied by “**CAUTION!**”, it means that you must follow an instruction to avoid damage to the product or other equipment.

Do not proceed until you clearly understand the instructions and all stated conditions are met.

Irritant substance


WEEE Compliance Statement

Valid for customers in EU countries

We are committed to being a good corporate citizen. As part of that commitment, we strive to maintain an environmentally conscious manufacturing operation. The European Union (EU) has enacted a directive on product recycling (Waste Electrical and Electronic Equipment, WEEE), and over time this directive will be implemented in the national laws of each EU member state.

Products falling under the scope of the WEEE Directive will be identified with a crossed over “wheeled bin” symbol on the product label, as indicated to the left. To forward the products to recycling or proper disposal, use an authorized collection system or return them to Pyrosequencing AB.

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Index

C

Computer Requirements	7
Contact Information	48

D

Delivered Items	
Analysis Software	7
Instrument	3
Vacuum Prep Workstation	9
Documentation	13

F

Function Test	
Instrument	12
Vacuum Prep Workstation	12

I

Install	
Analysis Software	7
Instrument	3
Vacuum Prep Workstation	9
Installation Test	
Instrument	12
Vacuum Prep Workstation	12
Instructions	13

L

Legal	1
-------	---

M

Move	
Instrument	3
Vacuum Prep Workstation	9

S

Safety	15
Label on the Instrument	18
Requirements	15
Safety in Other Languages	20
Sales Offices and Distributors	48
Service Agreements	14
Site Requirements	
Instrument	4
Vacuum Prep Workstation	10
Symbols	
Used in the User Documentation	18
Used on Product Labels	17

T

Technical Support	48
-------------------	----

U

Uninstall Analysis Software	8
Upgrade	
Analysis Software	7
Instrument Software	6
User Documentation	13

W

WEEE Compliance Statement	18
---------------------------	----

Contact Information

Pyrosequencing AB

Kungsgatan 76

SE-753 18 Uppsala

SWEDEN

Phone: +46 18 56 59 00

Fax: +46 18 59 19 22

E-mail: pyrosequencing@biotage.com

Website: www.biotagebio.com

Sales Offices and Distributors

Please visit www.biotagebio.com for contact information.

Technical Support

Europe: 1-pointsupport@eu.biotage.com

U.S. & Canada: 1-pointsupport@biotage.com

Rest of the world: Please contact your local distributor.

For phone numbers and more information, please visit www.biotagebio.com.

Revision History

Version	Revision date
AA	March 2008
AB	April 2008

April 2008

Code number 60-0276 AB