

## Primer Express Update

When importing a sequence into Primer Express® software v1.5a and v2.0, the software does not automatically convert U's into T's. To prevent the inadvertent deletion of the U bases, we recommend that the following steps be taken:

1. Copy the sequence into a word processing document (i.e. Microsoft Word)
2. In the Edit Menu select all, this will highlight all of your sequence.
3. Go back to the Edit Menu and Select "Replace"
4. The "Find and Replace" dialog box will open.
5. In the "Find what" text field type U.
6. In the "Replace with" text field type T.
7. Select replace all.
8. You can then copy and paste your sequence into Primer Express.