


Primer Express® P4 Enabler Overview


The Primer Express® P4 Enabler will allow users to run Primer Express® 2.0 on computers containing Pentium® 4 (Xeon) processors. This software updates two files in the original Primer Express® 2.0 release.

Installation Instructions


- 1) If installing from CD, copy the folder “Copy to Desktop” onto your computer’s desktop by dragging the folder icon. This is necessary, as the software will not run properly from the CD. Once copied, open the folder on your desktop by clicking on it. Proceed to step three.
- 2) If the software was downloaded from the Applied Biosystems web site, locate the PXP4.exe program. Proceed to step three.
- 3) Launch the software by clicking on the PXP4 icon. The following dialog box appears.




- 4) Click *OK*. The program will now begin searching for the installed file *primerexpress.exe*.



- 5) When the program has located *primerexpress.exe*, the following dialog box will be displayed.




- 6) Click **Yes** and the appropriate files will be updated. When the process has completed, exit the program by clicking on the **OK** button in the dialog shown below. The installation is complete.



Notes / Troubleshooting

Occasionally, the software is unable to find *primerexpress.exe*. If this is the case, the file will have to be located manually. The following dialog box will appear.



Manually locate the *primerexpress.exe* file and click *Open*. The dialog box shown in step 6 will then be displayed. Exit the program by clicking on the *OK* button. The installation is complete.

Legal Notice

This patch has been provided by PE Corporation (NY), through its Applied Bioystems Group – Informatics Division Customer Support Center (“Applied Biosystems”).

DISCLAIMER

This Supplemental Maintenance Release software (SMR) has been tested for specific functionality but has not been through a complete regression test. Due to the nature of SMRs, Applied Biosystems makes no representations regarding their use or performance and accepts no liability as a result of the use of this software. This software is provided "AS IS" and "WITH ALL FAULTS". Applied Biosystems provides no warranties either express or implied, including IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. In no event shall Applied Biosystems have any liability whatsoever for incidental and consequential damage as a result of the performance, use or operation of this software. The customer shall have sole responsibility for adequate protection and back-up of data used in conjunction with this SMR software.